# **EK-i/EW-i**

**Quick Start Guide – EK-i/EW-i Series  
Compact Balances**

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**AND**  
A&D Weighing

[andweighing.com](http://andweighing.com)

Quick Start Guide – EK-i/EW-i Series Compact Balances .....	1
Guide de démarrage rapide – Série EK-i/EW-i Balances compactes .....	11
Guía de inicio rápido – Serie EK-i/EW-i Básculas compactas .....	21

# EK-i/EW-i Series Compact Balances

## Quick Start Guide


 This is a hazard alert mark

---


## UNPACKING

---


When unpacking, check whether all the following items are included:

Main unit  
EW-150*i*  
EK-120*i* / EK-200*i* / EK-300*i*  
EK-410*i* / EK-610*i*


Main unit  
EW-1500*i* / EW-12K*i*  
EK-600*i* / EK-1200*i* / EK-2000*i*  
EK-3000*i* / EK-6000*i* / EK-12K*i*  
EK-4100*i* / EK-6100*i*

 Please confirm that the AC adapter is correct for your local voltage and receptacle type.

Quick Start Guide

---

# SETTING UP

---

### Setting up your balance


1. Place the weighing pan on the main unit as shown on the previous page.
2. Adjust the level of the balance using the leveling feet. Use the spirit level to confirm. The bubble should be in the center of the circle.
3. Calibrate your balance before use (see “CALIBRATION”).

### Balance location

To measure correctly, to keep the balance in good condition, and to prevent hazards; observe the following:

- Do not install the balance in locations that are subject to excessive dust, breezes, vibration, large temperature fluctuations, condensation, or that may have magnetic fields.
- Do not install the balance on a surface that is soft or that may cause the balance level to shift.
- Do not install the balance in direct sunshine.
- Do not install the balance near heaters or air conditioners.
- Do not use an unstable AC power source.
- Do not install the balance in a place where combustible or corrosive gases may exist.
- Allow the balance to reach equilibrium with the ambient temperature before use.
- Switch the power ON at least half an hour before use so that the balance can warm up.
- When the balance is installed for the first time, or the balance has been moved, carry out calibration as described in “CALIBRATION”.

**PART NAMES AND FUNCTIONS**

	<p>Turns the power on or off.</p>		<p>Outputs the weight data to a printer or computer, or enters a unit weight, 100% weight or other setting values to the balance.</p>
	<p>Hold down to enter the function setting mode. Weighing mode (EW-i only): Changes the weighing range (when rng 0 is selected). pcs mode: Enters the sample unit weight storing mode. % mode: Enters the 100% weight storing mode.</p>		<p>Switches the weighing unit (the weighing mode).</p>
			<p>Clears the display to zero.</p>


---

# CALIBRATION

---

### Calibration using a weight

1. Warm up the balance for at least half an hour with nothing on the pan.
2. Press and hold the calibration (CAL) switch until  appears, and release the switch.
3. The balance displays .

- ❖ To change the calibration weight value, proceed to step 4.
  - ❖ To use the calibration weight value in the balance memory, proceed to step 5.
4. Press the  key. The display shows the calibration weight value in “gram” that is stored in the balance. Use the following keys to change the value.

key To select the digit blinking to change.

key To set the value of the digit selected.

key To store the value and return to step 3.

key To cancel the value and return to step 3.

- Use a calibration weight of more than 2/3 of the capacity (of the highest range for EW-*i* series).
5. At step 3, pressing the  key weighs the zero-point value. Do not touch the pan during weighing.  
When the zero calibration is completed, the display shows the calibration weight value.
- To perform the zero calibration only and finish the procedure, press the  key.

6. Place the calibration weight with the same value as displayed on the pan. Press the  key to weigh it. Do not touch the pan during weighing.

7.  appears.

Remove the weight from the pan, and press the CAL switch or  key to return to the weighing mode.

#### Note:

**The value set in step 4 is stored in memory even after the power is switched off.**

**If the balance is to be moved to another place, set the gravity acceleration value of the area where the calibration using a weight is to be done, and calibrate the balance according to the procedure above. See the next section to set the value.**

---

## ERROR CODES

---

### Error codes

#### Overload error

e

Warning to indicate that an object beyond the balance capacity has been placed on the pan. Remove the object from the pan.

error5

Warning to indicate that the balance detects an unexpected load. Remove all objects from the pan or check that the weighing pan is correctly installed.

#### Range over notice

-e

This will be shown if the weight sensor receives strong force upward. Check if there is anything sandwiched around the weighing pan. There is a possibility that the weight sensor itself may have a failure.

#### CAL errors

Cal e

Warning to indicate that calibration has been canceled because the calibration weight is too heavy.

-Cal e

Warning to indicate that calibration has been canceled because the calibration weight is too light.

Check the weighing pan and the calibration weight. To return to the weighing mode, press the **MODE** key.

#### Low battery

lb0

Warning to indicate that the battery pack (OP-09) is discharged. Immediately stop using it, and recharge or use the AC adapter.

#### AC adapter error

Hb1

Warning to indicate that the output voltage of the AC adapter is too high. Check if the AC adapter is correct.

#### Stability error

error1

Warning to indicate that the weight value is not stable and the balance cannot display it. Prevent vibration and drafts. Press the **MODE** key to return to the weighing mode.

#### Internal errors (# = 2, 3, 4 or 6)

error#

The balance detects error state in the internal processing. With nothing on the weighing pan, turn the power off and on again. If the error persists, request service.

**If you cannot cancel an error or other errors have occurred, request service from the store where you purchased the balance or from your local A&D dealer.**

### Legal-for-Trade Model Information

#### NOTE FOR EC APPROVAL MODELS

EK-610i-EC, EK-6100i-EC, EW-150i-EC, EW-1500i-EC, EW-12Ki-EC

MODEL	EK-610i-EC	EK-6100i-EC
Class	II	II
Max. weight capacity	600g	6000g
Min. weight capacity	0.5g	5g
Min. display "d" g	0.01g	0.1g
Min. display "e" g	0.1g	1g
Max. weight capacity ct	3,000ct	30,000ct
Min. display "d" ct	0.05ct	0.5ct
Verification value "e" ct	0.5ct	5ct

MODEL	EW-150i-EC		EW-1500i-EC			EW-12Ki-EC		
Class	II		III			III		
Max. weight capacity	60g	150g	300g	600g	1500g	3kg	6kg	12kg
Min. weight capacity	0.4g	1g	2g	4g	10g	20g	40g	100g
Min. display "d" g	0.02g	0.05g	0.1g	0.2g	0.5g	1g	2g	5g
Min. display "e" g	0.02g	0.05g	0.1g	0.2g	0.5g	1g	2g	5g

\* In case of "Umin 0" (factory setting)

#### Note about function settings

Class	Item	Parameter	Description	
Func	rnG Range	0	Manual range change	Manual and Fixed ranges unavailable
		♦ 1	Automatic range change	
		2	Fixed to the lowest range	
		3	Fixed to the middle range	
		4	Fixed to the highest range	
	St-b Stability band width		Stability band width	Not available

**Note about Units:** When EU 2014/31/EU verification is performed only the following units will be available.

Unit Name	Unit Symbol	Unit Name	Unit Symbol
Grams	g	Pieces	pcs
Metric Carat	ct	Percent	%

## Quick Start Guide

### Note about Verification

#### Verification information for EU countries

The identification plate of the scale has a CE mark followed by the metrology mark (M + year two-digital; framed) and the number of the notified body. Such a scale can be started up and operated on-site provided that it is not connected to a peripheral device which is subject to legal control. Scales which are connected on site to a peripheral device which is subject to legal control must undergo a conformity assessment or verification. Only then, legal-for-trade use is permissible.

Scales without metrology mark or control marks may not be used in legal-for-trade applications. Scales with an accuracy – Class II, with 'n' greater than 3300 have been adjusted to the place of use and may not be used elsewhere.

As per legal requirements, the scale user must ensure that the scale is used for its intended purpose. This includes the observance of legal requirements, especially in case of changes, amendments and extensions.

#### Notes on verification

The verification of a scale is executed according to the respective statutory country regulations. The verification validity starts as soon as the device is put on the market.

---

## DISPOSAL OF THE DEVICE

---

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. To ensure that the device is disposed of correctly please recycle in accordance with local regulations. For further information contact the supplier from whom you purchased the product.

### EU Declaration of Conformity

This is a Declaration of Conformity for the instrument below.

<b>INSTRUMENT MODEL</b>	<i>EK-610i-EC &amp; EK-6100i-EC</i>
<b>NAME AND ADDRESS OF THE MANUFACTURER</b>	<b>A&amp;D Instruments Ltd</b> 24-26 Blacklands Way Abingdon Business Park Abingdon Oxfordshire OX14 1DY
<b>TYPE APPROVAL CERTIFICATE</b>	<b>T6969</b>
<b>NOTIFIED BODY ISSUING THE APPROVAL FOR MODULE B: EU TYPE EXAMINATION (Annex II .1 of 2014/31)</b>	<i>NMI 0122</i>
<b>NOTIFIED BODY ISSUING THE APPROVAL FOR MODULE D: CONFORMITY TO TYPE BASED ON QUALITY ASSURANCE OF THE PRODUCTION PROCESS (Annex II.2 of 2014/31)</b>	<i>SGS UK 0120</i>

Conformity to the following Directives is demonstrated by reference to the following harmonized standards or technical documentation:

<b>APPLICABLE DIRECTIVES</b>	<b>HARMONIZED STANDARDS OR OTHER TECHNICAL DOCUMENTATION</b>
Directive on non-automatic weighing instrument (2014/31/EU)	EN45501:2015
Directive on Electromagnetic Compatibility (2014/30/EU)	<i>EN 61326-1:2013</i> <i>Electrical equipment for measurement control and laboratory use – EMC requirements Part 1: General Requirements</i>
Directive on Electrical Equipment designed for use within certain voltage limits (2014/35/EU)	<i>EN 60950-1:2006 +A11:2009+A1:2010+A12:2011+A2:2013</i> <i>Safety of Information Technology Equipment</i>
Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS II) (2011/65/EU)	EN50581: 2012

This Declaration of Conformity is issued under the sole responsibility of the manufacturer:

Signed for and on behalf of:	<b>A&amp;D Instruments Ltd</b>
Place and date of issue:	Abingdon, UK, 1 <sup>st</sup> October 2018
Name, function:	J. Ghuman, Managing Director
Signature:	

## EU Declaration of Conformity

This is a Declaration of Conformity for the instrument below.

<b>INSTRUMENT MODEL</b>	<i>EW-150i-EC, EW-1500i-EC, EW-12Ki-EC</i>
<b>NAME AND ADDRESS OF THE MANUFACTURER</b>	<b>A&amp;D Instruments Ltd</b> <b>24-26 Blacklands Way</b> <b>Abingdon Business Park</b> <b>Abingdon</b> <b>Oxfordshire</b> <b>OX14 1DY</b>
<b>TYPE APPROVAL CERTIFICATE</b>	<b>T6391</b>
<b>NOTIFIED BODY ISSUING THE APPROVAL FOR MODULE B: EU TYPE EXAMINATION (Annex II .1 of 2014/31)</b>	<i>NMI 0122</i>
<b>NOTIFIED BODY ISSUING THE APPROVAL FOR MODULE D: CONFORMITY TO TYPE BASED ON QUALITY ASSURANCE OF THE PRODUCTION PROCESS (Annex II.2 of 2014/31)</b>	<i>SGS UK 0120</i>

Conformity to the following Directives is demonstrated by reference to the following harmonized standards or technical documentation

<b>APPLICABLE DIRECTIVES</b>	<b>HARMONIZED STANDARDS OR OTHER TECHNICAL DOCUMENTATION</b>
Directive on non-automatic weighing instrument (2014/31/EU)	EN45501:2015
Directive on Electromagnetic Compatibility (2014/30/EU)	<i>EN 61326-1:2013</i> <i>Electrical equipment for measurement control and laboratory use – EMC requirements Part 1: General Requirements</i>
Directive on Electrical Equipment designed for use within certain voltage limits (2014/35/EU)	<i>EN 60950-1:2006 +A11:2009+A1:2010+A12:2011+A2:2013</i> <i>Safety of Information Technology Equipment</i>
Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS II) (2011/65/EU)	EN50581: 2012

This Declaration of Conformity is issued under the sole responsibility of the manufacturer:

Signed for and on behalf of:	<b>A&amp;D Instruments Ltd</b>
Place and date of issue:	Abingdon, UK, 1 <sup>st</sup> October 2018
Name, function:	J. Ghuman, Managing Director
Signature:	

---

### EU Declaration of Conformity

---

A&D Instruments Ltd. hereby declare that the following Weighing product conforms to the requirements of the council directives on ...

**Electromagnetic Compatibility (EMC) 2014/30/EU,  
Low Voltage Equipment (LVD) 2014/35/EU and  
Restriction of the use of certain Hazardous Substances (RoHS) 2011/65/EU**

provided that they bear the CE mark of conformity.

**Model/Series....EK-i/EW-i Series**

Standards applicable:

**EN 61326-1:2013**

Electrical equipment for measurement, control and laboratory use -EMC requirements  
Part 1: General requirements

**EN 60950-1:2006/A11:2009/A1:2010/A12:2011/A2:2013**

Safety of Information Technology Equipment

**EN 50581:2012**

Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

This Declaration of Conformity is issued under the sole responsibility of the manufacturer:

A&D Instruments, 24-26 Blacklands Way, Abingdon, Oxfordshire, OX14 1DY, England

CE Mark first applied 08 July 2002

Signed for A&D Instruments on 1<sup>st</sup> October 2018

A handwritten signature in black ink, appearing to read 'J. Ghuman', is located below the text of the declaration.

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